

For Research Use Only

Mouse Il28a Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85590-3-PBS



Basic Information

Catalog Number:

85590-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001024673.2

GeneID (NCBI):

330496

UNIPROT ID:

Q4VK74

Full Name:

interleukin 28A

Calculated MW:

22 kDa

Purification Method:

Protein A purification

CloneNo.:

242247B2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01988-1, MOUSE IL-28A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85590-3-PBS. Detection antibody: 85590-2-PBS. Standard: Eg3169. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01988-2, MOUSE IL-28A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 98388-1-PBS. Detection antibody: 85590-3-PBS. Standard: Eg3169. Range: 0.313-20 ng/mL