

For Research Use Only

PMP2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85594-4-PBS



Basic Information

Catalog Number:

85594-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3411

GenBank Accession Number:

BC034997

GeneID (NCBI):

5375

UNIPROT ID:

P02689

Full Name:

peripheral myelin protein 2

Calculated MW:

132 aa, 15 kDa

Purification Method:

Protein A purification

CloneNo.:

243016E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01991-2, PMP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85594-4-PBS. Detection antibody: 85594-3-PBS. Standard: Ag3411. Range: 0.781-50 ng/mL