

For Research Use Only

PDIA5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85599-2-PBS



Basic Information

Catalog Number:

85599-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG31479

GenBank Accession Number:

BC001625

GeneID (NCBI):

10954

UNIPROT ID:

Q14554

Full Name:

protein disulfide isomerase family A, member 5

Calculated MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

242971E6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01965-2, PDIR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85599-3-PBS. Detection antibody: 85599-2-PBS. Standard: Ag31479. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP01965-1, PDIR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85599-1-PBS. Detection antibody: 85599-2-PBS. Standard: Ag31479. Range: 1.563-100 ng/mL.