

For Research Use Only

# JAM2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85601-1-PBS



## Basic Information

Catalog Number:

85601-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_021219.3

GeneID (NCBI):

58494

UNIPROT ID:

P57087-1

Full Name:

junctional adhesion molecule 2

Calculated MW:

33 kDa

Purification Method:

Protein A purification

CloneNo.:

242617G8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01977-1, JAM2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85601-2-PBS. Detection antibody: 85601-1-PBS. Standard: Eg2914. Range: 0.313-40 ng/mL