

For Research Use Only

GLUD2 Recombinant antibody, PBS Only

Catalog Number: 85602-5-PBS



Basic Information

Catalog Number:

85602-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5748

GenBank Accession Number:

BC050732

GeneID (NCBI):

2747

UNIPROT ID:

P49448

Full Name:

glutamate dehydrogenase 2

Calculated MW:

61 kDa

Observed MW:

50-60 kDa

Purification Method:

Protein A purification

CloneNo.:

243041D11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85602-5-RR (GLUD2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85602-5-PBS in a different storage buffer formulation.