

For Research Use Only

GLUD2 Recombinant antibody

Catalog Number: 85602-5-RR



Basic Information

Catalog Number:

85602-5-RR

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5748

GenBank Accession Number:

BC050732

GeneID (NCBI):

2747

UNIPROT ID:

P49448

Full Name:

glutamate dehydrogenase 2

Calculated MW:

61 kDa

Observed MW:

50-60 kDa

Purification Method:

Protein A purification

CloneNo.:

243041D11

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HeLa cells, HepG2 cells, A431 cells

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85602-5-RR (GLUD2 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.