

For Research Use Only

GLUD2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85602-6-PBS



Basic Information

Catalog Number:

85602-6-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5748

GenBank Accession Number:

BC050732

GeneID (NCBI):

2747

UNIPROT ID:

P49448

Full Name:

glutamate dehydrogenase 2

Calculated MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

243041G12

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01985-3, Human GLUD2 Recombinant Matched Antibody Pair - PBS only. 85602-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag5748. 85602-1-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL.

The mean GLUD2 concentration was determined to be 21.35 ng/mL in HeLa cell extract based on a 1.80 mg/mL extract load.