

For Research Use Only

# EIF3I Recombinant antibody, PBS Only (Capture)

Catalog Number: 85603-3-PBS



## Basic Information

Catalog Number:

85603-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1847

GenBank Accession Number:

BC000413

GeneID (NCBI):

8668

UNIPROT ID:

Q13347

Full Name:

eukaryotic translation initiation factor 3, subunit I

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

243040B1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01969-2, EIF3I Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85603-3-PBS. Detection antibody: 85603-2-PBS. Standard: Ag1847. Range: 0.781-100 ng/mL.