

For Research Use Only

CHGA Recombinant antibody, PBS Only (Detector)

Catalog Number: 85605-1-PBS



Basic Information

Catalog Number:

85605-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19780

GenBank Accession Number:

BC006459

GeneID (NCBI):

1113

UNIPROT ID:

P10645

Full Name:

chromogranin A (parathyroid secretory protein 1)

Calculated MW:

51 kDa

Purification Method:

Protein A purification

CloneNo.:

243048E8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01975-2, Chromogranin A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85605-2-PBS. Detection antibody: 85605-1-PBS. Standard: Ag19780. Range: 3.125-100 ng/mL.