

For Research Use Only

# NUP133 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85609-2-PBS



## Basic Information

Catalog Number:

85609-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3100

GenBank Accession Number:

BC020107

GeneID (NCBI):

55746

UNIPROT ID:

Q8WUM0

Full Name:

nucleoporin 133kDa

Calculated MW:

1156 aa, 129 kDa

Purification Method:

Protein A purification

CloneNo.:

243092G3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01982-1, NUP133 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85609-2-PBS. Detection antibody: 85609-3-PBS. Standard: Ag3100. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01982-2, NUP133 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85609-2-PBS. Detection antibody: 85609-1-PBS. Standard: Ag3100. Range: 0.781-100 ng/mL