

For Research Use Only

PRNP Recombinant antibody, PBS Only (Detector)

Catalog Number: 85611-2-PBS



Basic Information

Catalog Number:

85611-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC012844

GeneID (NCBI):

5621

UNIPROT ID:

F7VJQ1

Full Name:

prion protein

Purification Method:

Protein A purification

CloneNo.:

243104B12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01974-2, PrP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85611-3-PBS. Detection antibody: 85611-2-PBS. Standard: Eg3359. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01974-1, PrP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85611-1-PBS. Detection antibody: 85611-2-PBS. Standard: Eg3359. Range: 0.781-100 ng/mL.