

For Research Use Only

# Mouse Tim4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85616-2-PBS



## Basic Information

Catalog Number:

85616-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_178759

GeneID (NCBI):

276891

UNIPROT ID:

Q6U7R4

Full Name:

T-cell immunoglobulin and mucin domain containing 4

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

242616A7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01983-1, MOUSE TIMD4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85616-2-PBS. Detection antibody: 85616-1-PBS. Standard: Eg3307. Range: 0.781-100 ng/mL.