

For Research Use Only

MEOX1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85618-3-PBS



Basic Information

Catalog Number:

85618-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12910

GenBank Accession Number:

BC069474

GeneID (NCBI):

4222

UNIPROT ID:

P50221

Full Name:

mesenchyme homeobox 1

Calculated MW:

254 aa, 28 kDa

Purification Method:

Protein A purification

CloneNo.:

242909A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02006-2, MEOX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85618-3-PBS. Detection antibody: 85618-2-PBS. Standard: Ag12910. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP02006-1, MEOX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85618-1-PBS. Detection antibody: 85618-3-PBS. Standard: Ag12910. Range: 0.625-80 ng/mL