

For Research Use Only

# PLEK2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85619-3-PBS



## Basic Information

Catalog Number:

85619-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2297

GenBank Accession Number:

BC001226

GeneID (NCBI):

26499

UNIPROT ID:

Q9NYT0

Full Name:

pleckstrin 2

Calculated MW:

353 aa, 40 kDa

Purification Method:

Protein A purification

CloneNo.:

242987E5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP02002-1, PLEK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85619-3-PBS. Detection antibody: 85619-2-PBS. Standard: Ag2297. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP02002-2, PLEK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85619-3-PBS. Detection antibody: 85619-1-PBS. Standard: Ag2297. Range: 0.781-100 ng/mL.