

For Research Use Only

Mouse Ngfr Recombinant antibody, PBS Only (Capture)

Catalog Number: 85622-3-PBS



Basic Information

Catalog Number:

85622-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_033217.3

GeneID (NCBI):

18053

UNIPROT ID:

Q9Z0W1-1

Full Name:

nerve growth factor receptor (TNFR superfamily, member 16)

Calculated MW:

46 kDa

Purification Method:

Protein A purification

CloneNo.:

242634F11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02004-1, MOUSE NGFR/TNFRSF16 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85622-3-PBS. Detection antibody: 85622-2-PBS. Standard: Eg2651. Range: 0.781-50 ng/mL