

For Research Use Only

Mouse Tslp Recombinant antibody, PBS Only (Detector)

Catalog Number:85623-2-PBS



Basic Information

Catalog Number:

85623-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_021367.2

GeneID (NCBI):

53603

UNIPROT ID:

Q9JIE6

Full Name:

thymic stromal Lymphopoietin

Calculated MW:

16KD

Purification Method:

Protein A purification

CloneNo.:

242914A9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01963-2, MOUSE Tslp Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85623-3-PBS. Detection antibody: 85623-2-PBS. Standard: Eg3312. Range: 0.625-40 ng/mL.