

For Research Use Only

Mouse Cd300lg Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number: 85624-3-PBS



Basic Information

Catalog Number: 85624-3-PBS	GenBank Accession Number: NM_001160712.1	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 52685	CloneNo.: 242032A9
Isotype: IgG	UNIPROT ID: NP_001154183.1	
Immunogen Catalog Number: EG2731	Full Name: CD300 antigen like family member G	
	Calculated MW: 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01979-1, MOUSE Cd300lg Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85450-1-PBS. Detection antibody: 85450-2-PBS. Standard: Eg2731. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01979-2, MOUSE Cd300lg Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85450-1-PBS. Detection antibody: 85450-3-PBS. Standard: Eg2731. Range: 0.625-80 ng/mL