

For Research Use Only

SOCS4 Recombinant antibody, PBS Only

Catalog Number: 85626-1-PBS



Basic Information

Catalog Number:

85626-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35370

GenBank Accession Number:

BC060790

GeneID (NCBI):

122809

UNIPROT ID:

Q8WXH5

Full Name:

suppressor of cytokine signaling 4

Calculated MW:

51 kDa

Observed MW:

51 kDa

Purification Method:

Protein A purification

CloneNo.:

243060E3

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

SOCS4 contains a SH2 domain and a SOCS BOX domain, and belongs to Suppressors of cytokine signaling (SOCS) family. SOCS family members are known to be cytokine-inducible negative regulators of cytokine signaling.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85626-1-RR (SOCS4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85626-1-PBS in a different storage buffer formulation.