

For Research Use Only

# SOCS4 Recombinant antibody

Catalog Number: 85626-1-RR



## Basic Information

Catalog Number:

85626-1-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35370

GenBank Accession Number:

BC060790

GeneID (NCBI):

122809

UNIPROT ID:

Q8WXH5

Full Name:

suppressor of cytokine signaling 4

Calculated MW:

51 kDa

Observed MW:

51 kDa

Purification Method:

Protein A purification

CloneNo.:

243060E3

Recommended Dilutions:

WB 1:5000-1:20000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : MCF-7 cells, Jurkat cells, U-87 MG cells, HeLa cells, HEK-293T cells

## Background Information

SOCS4 contains a SH2 domain and a SOCS BOX domain, and belongs to Suppressors of cytokine signaling (SOCS) family. SOCS family members are known to be cytokine-inducible negative regulators of cytokine signaling.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85626-1-RR (SOCS4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.