

For Research Use Only

# HERPUD1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85627-2-PBS



## Basic Information

Catalog Number:

85627-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34891

GenBank Accession Number:

BC008320

GeneID (NCBI):

9709

UNIPROT ID:

Q15011

Full Name:

homocysteine-inducible,  
endoplasmic reticulum stress-  
inducible, ubiquitin-like domain  
member 1

Calculated MW:

44 kDa

Purification Method:

Protein A purification

CloneNo.:

243054A5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02032-2, HERPUD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85627-3-PBS. Detection antibody: 85627-2-PBS. Standard: Ag34891. Range: 0.781-100 ng/mL.