

For Research Use Only

HERPUD1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number:85627-3-PBS



Basic Information

Catalog Number:

85627-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34891

GenBank Accession Number:

BC008320

GeneID (NCBI):

9709

UNIPROT ID:

Q15011

Full Name:

homocysteine-inducible,
endoplasmic reticulum stress-
inducible, ubiquitin-like domain
member 1

Calculated MW:

44 kDa

Purification Method:

Protein A purification

CloneNo.:

24305489

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MPO2032-1, HERPUD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85627-4-PBS. Detection antibody: 85627-3-PBS. Standard: Ag34891. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MPO2032-2, HERPUD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85627-3-PBS. Detection antibody: 85627-2-PBS. Standard: Ag34891. Range: 0.781-100 ng/mL