

For Research Use Only

# Rat Lumican Recombinant antibody, PBS Only (Capture)

Catalog Number: 85629-4-PBS



## Basic Information

Catalog Number:

85629-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3410

GenBank Accession Number:

NM\_031050.1

GeneID (NCBI):

81682

UNIPROT ID:

P51886

Full Name:

lumican

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

242298D6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02015-2, RAT Lumican Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85629-4-PBS. Detection antibody: 85629-3-PBS. Standard: Eg3410. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP02015-1, RAT Lumican Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85629-4-PBS. Detection antibody: 85629-2-PBS. Standard: Eg3410. Range: 0.781-100 ng/mL