

For Research Use Only

Rat Osteoactivin/GPNMB Recombinant antibody, PBS Only (Capture)

Catalog Number: 85632-4-PBS



Basic Information

Catalog Number:

85632-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3283

GenBank Accession Number:

NM_133298.1

GeneID (NCBI):

113955

UNIPROT ID:

Q9QZF6

Full Name:

glycoprotein (transmembrane) nmb

Calculated MW:

64 kDa

Purification Method:

Protein A purification

CloneNo.:

242601A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02036-2, RAT GPNMB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85632-4-PBS. Detection antibody: 85632-3-PBS. Standard: Eg3283. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP02036-1, RAT GPNMB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85632-4-PBS. Detection antibody: 85632-2-PBS. Standard: Eg3283. Range: 1.563-200 ng/mL.