

For Research Use Only

# RPA1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85633-3-PBS



## Basic Information

Catalog Number:

85633-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14543

GenBank Accession Number:

BC017115

GeneID (NCBI):

25885

UNIPROT ID:

O95602

Full Name:

polymerase (RNA) I polypeptide A,  
194kDa

Calculated MW:

1720 aa, 195 kDa

Purification Method:

Protein A purification

CloneNo.:

243082F5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP02120-2, RPA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85633-3-PBS. Detection antibody: 85633-4-PBS. Standard: Ag14543. Range: 0.781-100 ng/mL.