

For Research Use Only

RPA1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85633-5-PBS



Basic Information

Catalog Number:

85633-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14543

GenBank Accession Number:

BC017115

GeneID (NCBI):

25885

UNIPROT ID:

O95602

Full Name:

polymerase (RNA) I polypeptide A,
194kDa

Calculated MW:

1720 aa, 195 kDa

Purification Method:

Protein A purification

CloneNo.:

243082B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02120-1, RPA1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85633-5-PBS. Detection antibody: 85633-4-PBS. Standard: Ag14543. Range: 0.781-25 ng/mL