

For Research Use Only

CDO1 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85636-2-PBS



Basic Information

Catalog Number:

85636-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3294

GenBank Accession Number:

BC024241

GeneID (NCBI):

1036

UNIPROT ID:

Q16878

Full Name:

cysteine dioxygenase, type I

Calculated MW:

200 aa, 23 kDa

Purification Method:

Protein A purification

CloneNo.:

250039D9

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02053-1, CDO1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85636-3-PBS. Detection antibody: 85636-2-PBS. Standard: Ag3294. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP02053-2, CDO1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85636-2-PBS. Detection antibody: 85636-1-PBS. Standard: Ag3294. Range: 0.781-100 ng/mL.