

For Research Use Only

FARP1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85638-2-PBS



Basic Information

Catalog Number:

85638-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21225

GenBank Accession Number:

BC071592

GeneID (NCBI):

10160

UNIPROT ID:

Q9Y4F1

Full Name:

FERM, RhoGEF (ARHGEF) and pleckstrin domain protein 1 (chondrocyte-derived)

Calculated MW:

1076 aa, 122 kDa

Purification Method:

Protein A purification

CloneNo.:

243132B6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02060-1, FARP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85638-3-PBS. Detection antibody: 85638-2-PBS. Standard: Ag21225. Range: 0.469-60 ng/mL



Cytometric bead array standard curve of MP02060-2, FARP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85638-1-PBS. Detection antibody: 85638-2-PBS. Standard: Ag21225. Range: 0.469-60 ng/mL