

For Research Use Only

# PFKFB2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85639-5-PBS



## Basic Information

Catalog Number:

85639-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG32347

GenBank Accession Number:

BC075076

GeneID (NCBI):

5208

UNIPROT ID:

O60825

Full Name:

6-phosphofructo-2-kinase/fructose-2,6-biphosphatase 2

Calculated MW:

505 aa, 58 kDa

Purification Method:

Protein A purification

CloneNo.:

243136A8

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP02061-3, Human PFKFB2 Recombinant Matched Antibody Pair - PBS only. 85639-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag32347. 85639-5-PBS was HRP conjugated as the detection antibody. Range: 0.313-20 ng/mL



The mean PFKFB2 concentration was determined to be 22.84 ng/mL in HepG2 cell extract based on a 1.8 mg/mL extract load and 24.89 ng/mL in HeLa cell extract based on a 1.8 mg/mL extract load.