

For Research Use Only

# SF3B1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85640-2-PBS



## Basic Information

Catalog Number:

85640-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26689

GenBank Accession Number:

NM\_012433

GeneID (NCBI):

23451

UNIPROT ID:

O75533

Full Name:

splicing factor 3b, subunit 1, 155kDa

Calculated MW:

146 kDa

Purification Method:

Protein A purification

CloneNo.:

243145B4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02058-2, SF3B1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85640-3-PBS. Detection antibody: 85640-2-PBS. Standard: Ag26689. Range: 0.781-100 ng/mL.