

For Research Use Only

THOC1 Recombinant antibody, PBS Only

Catalog Number: 85642-5-PBS



Basic Information

Catalog Number: 85642-5-PBS	GenBank Accession Number: BC010381	Purification Method: Protein A purification
Concentration: 1000 µg/ml	GeneID (NCBI): 9984	CloneNo.: 243161E9
Source: Rabbit	UNIPROT ID: Q96FV9	
Isotype: IgG	Full Name: THO complex 1	
Immunogen Catalog Number: AG1352	Calculated MW: 84 kDa	
	Observed MW: 76 kDa	

Applications

Tested Applications:
WB, IF/ICC, Indirect ELISA

Species Specificity:
human, mouse, rat

Background Information

THOC1 (THO complex subunit 1), also known as Tho1, P84, HPR1 or P84N5, is a 657 amino acid nuclear matrix protein and is evolutionarily conserved from yeast to humans. THOC1 contains one death domain and is a component of the heteromultimeric THO/TREX (transcription/export) complex along with THOC2, THOC3, BAT1 and ALY. The THO/TREX complex is recruited to transcribed genes and travels along with RNA polymerase II (Pol II) during elongation, coupling elongating Pol II with RNA splicing and export factors. THOC1 is expressed at high levels in breast cancer cells and at relatively low levels in normal epithelia. A reduction of THOC1 in cancer cell lines results in promoted cancer cell apoptosis and reduced cell proliferation. This suggests that cancer cells are dependent on the high levels of THOC1 expression and therefore THOC1 may be a good target for cancer therapy.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 85642-5-RR (THOC1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85642-5-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using THOC1 antibody (85642-5-RR, Clone: 243161E9) at dilution of 1:800 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 85642-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 85642-5-RR against Human THOC1 were performed. The affinity constant is below 1 pM.