

For Research Use Only

CXCL4/PF4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85647-3-PBS



Basic Information

Catalog Number:

85647-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002619.4

GeneID (NCBI):

5196

UNIPROT ID:

P02776

Full Name:

platelet factor 4

Calculated MW:

11 kDa

Purification Method:

Protein A purification

CloneNo.:

243087E3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02040-2, CXCL4/PF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85647-3-PBS. Detection antibody: 98410-2-PBS. Standard: Eg2129. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP02040-1, CXCL4/PF4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85647-3-PBS. Detection antibody: 85647-2-PBS. Standard: Eg2129. Range: 0.781-100 ng/mL