

For Research Use Only

SNRPA1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85649-5-PBS



Basic Information

Catalog Number:

85649-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10790

GenBank Accession Number:

BC022816

GeneID (NCBI):

6627

UNIPROT ID:

P09661

Full Name:

small nuclear ribonucleoprotein

polypeptide A'

Calculated MW:

255 aa, 28 kDa

Purification Method:

Protein A purification

CloneNo.:

243097E9

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02027-3, Human SNRPA1 Recombinant Matched Antibody Pair - PBS only. 85649-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag10790. 85649-5-PBS was HRP conjugated as the detection antibody. Range: 39.1-2500 pg/mL.



The mean SNRPA1 concentration was determined to be 211.7 ng/mL in HeLa cell extract based on a 2.1 mg/mL extract load and 334.5 ng/mL in Jurkat cell extract based on a 1.9 mg/mL extract load.