

For Research Use Only

IQGAP1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85651-4-PBS



Basic Information

Catalog Number:

85651-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG19321

GenBank Accession Number:

BC151834

GeneID (NCBI):

8826

ENSEMBL Gene ID:

ENSG00000140575

UNIPROT ID:

P46940

Full Name:

IQ motif containing GTPase activating protein 1

Calculated MW:

189 kDa

Purification Method:

Protein A purification

CloneNo.:

242206F10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02009-1, IQGAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85651-5-PBS. Detection antibody: 85651-4-PBS. Standard: Ag19321. Range: 1.563-200 ng/mL



Cytometric bead array standard curve of MP02009-2, IQGAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85651-3-PBS. Detection antibody: 85651-4-PBS. Standard: Ag19321. Range: 1.563-200 ng/mL