

For Research Use Only

# BZW1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85655-3-PBS



## Basic Information

Catalog Number:

85655-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13830

GenBank Accession Number:

BC001804

GeneID (NCBI):

9689

UNIPROT ID:

Q7L1Q6

Full Name:

basic leucine zipper and W2 domains

1

Calculated MW:

353 aa, 41 kDa

Purification Method:

Protein A purification

CloneNo.:

243162F6

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP02022-2, BZW1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85655-3-PBS. Detection antibody: 85655-2-PBS. Standard: Ag13830. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP02022-1, BZW1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85655-3-PBS. Detection antibody: 85655-1-PBS. Standard: Ag13830. Range: 0.781-100 ng/mL.