

For Research Use Only

FLT3LG Recombinant antibody, PBS Only (Detector)

Catalog Number: 85656-3-PBS



Basic Information

Catalog Number:

85656-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001459.3

GeneID (NCBI):

2323

UNIPROT ID:

P49771-1

Full Name:

fms-related tyrosine kinase 3 ligand

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

242884D3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02023-1, FLT3 Ligand Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85656-4-PBS. Detection antibody: 85656-3-PBS. Standard: Eg2926. Range: 0.78-100 ng/mL