

For Research Use Only

Mouse BAFFR/TNFRSF13C Recombinant antibody, PBS Only (Capture)

Catalog Number: 85658-1-PBS



Basic Information

Catalog Number:

85658-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_028075.2

GeneID (NCBI):

72049

UNIPROT ID:

Q9D8D0-1

Full Name:

tumor necrosis factor receptor
superfamily, member 13c

Calculated MW:

19kDa

Purification Method:

Protein A purification

CloneNo.:

243120F4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02045-2, MOUSE TNFRSF13C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85658-1-PBS. Detection antibody: 85658-2-PBS. Standard: Eg1840. Range: 0.781-100 ng/mL.