

For Research Use Only

# Mouse Complement Component C5a Recombinant antibody, PBS Only (Detector)

Catalog Number: 85659-1-PBS



## Basic Information

Catalog Number:

85659-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_010406.2

GeneID (NCBI):

15139

UNIPROT ID:

P06684

Full Name:

hemolytic complement

Calculated MW:

189kDa

Purification Method:

Protein A purification

CloneNo.:

242757C10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02047-1, MOUSE C5a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85659-2-PBS. Detection antibody: 85659-1-PBS. Standard: Eg3332. Range: 0.781-100 ng/mL.