

For Research Use Only

# Mouse Il2rb (CD122) Recombinant antibody, PBS Only (Capture)

Catalog Number: 85660-3-PBS



## Basic Information

Catalog Number:

85660-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_008368.4

GeneID (NCBI):

16185

UNIPROT ID:

P16297

Full Name:

interleukin 2 receptor, beta chain

Calculated MW:

61 kDa

Purification Method:

Protein A purification

CloneNo.:

243118G2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02046-2, MOUSE CD122/IL-2RB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85660-3-PBS. Detection antibody: 85660-2-PBS. Standard: Eg2709. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP02046-1, MOUSE CD122/IL-2RB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85660-3-PBS. Detection antibody: 85660-1-PBS. Standard: Eg2709. Range: 0.781-100 ng/mL