

For Research Use Only

PLOD1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85662-3-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 85662-3-PBS | GenBank Accession Number: BC016657 | Purification Method: Protein A purification |
| Concentration: 1 mg/ml | GeneID (NCBI): 5351 | CloneNo.: 242403F8 |
| Source: Rabbit | UNIPROT ID: Q02809 | |
| Isotype: IgG | Full Name: procollagen-lysine 1, 2-oxoglutarate 5-dioxygenase 1 | |
| Immunogen Catalog Number: AG31626 | Calculated MW: 727 aa, 84 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02019-2, PLOD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85662-4-PBS. Detection antibody: 85662-3-PBS. Standard: Ag31626. Range: 0.781-100 ng/mL.