

For Research Use Only

RIOK1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85668-2-PBS



Basic Information

Catalog Number:

85668-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11057

GenBank Accession Number:

BC006104

GeneID (NCBI):

83732

UNIPROT ID:

Q9BRS2

Full Name:

RIO kinase 1 (yeast)

Calculated MW:

568 aa, 66 kDa

Purification Method:

Protein A purification

CloneNo.:

243091E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02010-2, R1OK1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85668-2-PBS. Detection antibody: 85668-1-PBS. Standard: Ag11057. Range: 0.781-100 ng/mL.