

For Research Use Only

SKP2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85670-2-PBS



Basic Information

Catalog Number:

85670-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6826

GenBank Accession Number:

BC001441

GeneID (NCBI):

6502

UNIPROT ID:

Q13309

Full Name:

S-phase kinase-associated protein 2 (p45)

Calculated MW:

48 kDa

Purification Method:

Protein A purification

CloneNo.:

243080G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02007-2, SKP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85670-3-PBS. Detection antibody: 85670-2-PBS. Standard: Ag6826. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP02007-1, SKP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85670-1-PBS. Detection antibody: 85670-2-PBS. Standard: Ag6826. Range: 0.781-100 ng/mL.