

For Research Use Only

Mouse Il15ra Recombinant antibody, PBS Only (Capture)

Catalog Number: 85671-4-PBS



Basic Information

Catalog Number:

85671-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_008358.2

GeneID (NCBI):

16169

UNIPROT ID:

Q60819-1

Full Name:

interleukin 15 receptor, alpha chain

Calculated MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

242618G11

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02016-3, Mouse IL-15RA Recombinant Matched Antibody Pair - PBS only. 85671-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2687. 85671-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL

The mean IL-15RA concentration was determined to be 1,332.06 ng/mL in mouse liver tissue extract based on a 3.7 mg/mL extract load.