

For Research Use Only

Rat IL-1RA/IL-1RN Recombinant antibody, PBS Only (Capture)

Catalog Number: 85672-1-PBS



Basic Information

Catalog Number:	85672-1-PBS	GenBank Accession Number:	NM_022194.2	Purification Method:	Protein A purification
Concentration:	1 mg/ml	GeneID (NCBI):	60582	CloneNo.:	242317E10
Source:	Rabbit	UNIPROT ID:	P25086		
Isotype:	IgG	Full Name:	interleukin 1 receptor antagonist		
Immunogen Catalog Number:	EG3206	Calculated MW:	20 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
rat

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02044-1, RAT IL-1RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85672-1-PBS. Detection antibody: 85672-3-PBS. Standard: Eg3206. Range: 1.563-200 ng/mL.