

For Research Use Only

# Rat IL-1RA/IL-1RN Recombinant antibody, PBS Only (Detector)

Catalog Number: 85672-2-PBS



## Basic Information

Catalog Number:

85672-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3206

GenBank Accession Number:

NM\_022194.2

GeneID (NCBI):

60582

UNIPROT ID:

P25086

Full Name:

interleukin 1 receptor antagonist

Calculated MW:

20 kDa

Purification Method:

Protein A purification

CloneNo.:

242317A10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02044-2, RAT IL-1RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85672-3-PBS. Detection antibody: 85672-2-PBS. Standard: Eg3206. Range: 1.563-200 ng/mL.