

For Research Use Only

# Rat Osteopontin Recombinant antibody, PBS Only (Detector)

Catalog Number: 85673-2-PBS



## Basic Information

Catalog Number:

85673-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3247

GenBank Accession Number:

NM\_012881.2

GeneID (NCBI):

25353

UNIPROT ID:

P08721

Full Name:

secreted phosphoprotein 1

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

24265483

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02043-1, RAT Osteopontin/OPN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85673-3-PBS. Detection antibody: 85673-2-PBS. Standard: Eg3247. Range: 0.156-20 ng/mL