

For Research Use Only

# CXCL16 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85680-2-PBS



## Basic Information

Catalog Number:

85680-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC017588

GeneID (NCBI):

58191

UNIPROT ID:

Q9H2A7

Full Name:

chemokine (C-X-C motif) ligand 16

Calculated MW:

273 aa, 30 kDa

Purification Method:

Protein A purification

CloneNo.:

242984H2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP02039-1, CXCL16 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85680-3-PBS. Detection antibody: 85680-2-PBS. Standard: Eg2165. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP02039-2, CXCL16 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85680-1-PBS. Detection antibody: 85680-2-PBS. Standard: Eg2165. Range: 0.781-100 ng/mL.