

For Research Use Only

CXCL16 Recombinant antibody, PBS Only

Catalog Number: 85680-4-PBS



Basic Information

Catalog Number: 85680-4-PBS	GenBank Accession Number: BC017588	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 58191	CloneNo.: 242984H8
Source: Rabbit	UNIPROT ID: Q9H2A7	
Isotype: IgG	Full Name: chemokine (C-X-C motif) ligand 16	
	Calculated MW: 273 aa, 30 kDa	
	Observed MW: 28-30 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human

Background Information

CXCL16 is a recently discovered cytokine belonging to the CXC chemokine family, which is synthesised in plasmacytoid dendritic cell as a transmembrane molecule. It exists in a transmembrane and soluble form. The transmembrane form of CXCL16 functions as an adhesion molecule for CXCR6-expressing cells, whereas the soluble form of CXCL16 mediates infiltration of circulating cells into sites of injury. CXCL16, has been proposed as an important pathogenic mediator in inflammatory diseases, including rheumatoid arthritis, glomerulonephritis, or prostate cancer. CXCL16 has been implicated in some forms of renal disease such as lupus nephritis and antglomerular basement membrane nephritis. CXCL16 also plays a pivotal role in the pathogenesis of angiotensin II-induced renal injury and fibrosis through regulation of macrophage and T cell infiltration and bone marrow-derived fibroblast accumulation.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



human testis tissue were subjected to SDS PAGE followed by western blot with 85680-4-RR (CXCL16 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 85680-4-PBS in a different storage buffer formulation.