

For Research Use Only

KTN1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85681-3-PBS



Basic Information

Catalog Number:

85681-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13854

GenBank Accession Number:

BC117132

GeneID (NCBI):

3895

UNIPROT ID:

Q86UP2

Full Name:

kinectin 1 (kinesin receptor)

Calculated MW:

1357 aa, 156 kDa

Purification Method:

Protein A purification

CloneNo.:

243067D10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02037-1, KTN1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85681-3-PBS. Detection antibody: 85681-2-PBS. Standard: Ag13854. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP02037-2, KTN1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85681-3-PBS. Detection antibody: 85681-1-PBS. Standard: Ag13854. Range: 3.125-100 ng/mL.