

For Research Use Only

# SERPINE1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85683-3-PBS



## Basic Information

Catalog Number:

85683-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC010860

GeneID (NCBI):

5054

UNIPROT ID:

P05121

Full Name:

serpin peptidase inhibitor, clade E  
(nexin, plasminogen activator  
inhibitor type 1), member 1

Calculated MW:

402 aa, 45 kDa

Purification Method:

Protein A purification

CloneNo.:

243115F8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP02035-2, Serpin E1/PAI-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85683-3-PBS. Detection antibody: 85683-2-PBS. Standard: Eg3798. Range: 1.563-200 ng/mL.