

For Research Use Only

SPINT2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85684-4-PBS



Basic Information

Catalog Number:

85684-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_021102.3

GeneID (NCBI):

10653

UNIPROT ID:

O43291-1

Full Name:

serine peptidase inhibitor, Kunitz type, 2

Calculated MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

243094F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02034-1, SPINT2/HAI-2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85684-4-PBS. Detection antibody: 85684-3-PBS. Standard: Eg3130. Range: 0.781-100 ng/mL.